

Saule<sup>TM</sup>

Rx Pre-dilatation PTCA Balloon Catheter

Specification

|                               |                                                  |
|-------------------------------|--------------------------------------------------|
| Balloon Diameter(mm)          | 1.5,1.75,2.0,2.25,2.5,2.75,3.0,3.25,3.5,3.75,4.0 |
| Balloon Length(mm)            | 9,10,12,14,15,16,20,24,25,30                     |
| Nominal Pressure(atm)         | 6                                                |
| Rated Burst Pressure(atm)     | 16                                               |
| Guidewire Compatibility(inch) | 0.014                                            |
| Catheter Length(mm)           | 1400                                             |
| Distal Outer Diameter(F)      | 2.7                                              |
| Lesion Entry Profile(F)       | 0.0165                                           |



Princeton Medical Scientific Co., Ltd  
Rm 503, Building C, No.8 Fukang West Street,  
Yushu Industrial Park, Guangzhou Science City,  
Guangzhou, 510663, China  
T: (86)20-82112462  
F: (86)20-82119811  
E: pulin@pulinmed.com  
www.pulinmed.com



Saule<sup>TM</sup>

Rx Pre-dilatation PTCA Balloon Catheter



Designed for all the lesions

www.pulinmed.com

## Application

1. Intended for balloon dilatation of the stenotic lesion of a coronary artery for the purpose of improving blood flow.
2. Used for stent post-dilatation in the coronary vasculature.

## Features

### Seamless bonding technology

Tip and guidewire seamless bonding enables excellent tracking and avoid fish mouth phenomenon

### Flexible and soft tip design


Quicker and easier to pass the severe stenotic and calcified lesion with the most minimized trauma

### Precision laser bonding technology


Stepless connection of the balloon and catheter provides superior pushability

### Multi- folded technology

Durable hydrophilic coating and lower entry profile provide successful crossability for complex lesion



## Before and after PCI procedure



Before



After





### Specification

| Shaft Length 1400mm, Guide Wire 0.014", RX |                    |                  |                  |                  |                  |                  |
|--------------------------------------------|--------------------|------------------|------------------|------------------|------------------|------------------|
| Balloon Diameter (mm)                      | Balloon Length(mm) |                  |                  |                  |                  |                  |
|                                            | 9                  | 10               | 12               | 14               | 15               | 16               |
| 1.5                                        | 04 015 009 140 2   | 04 015 010 140 2 | 04 015 012 140 2 | 04 015 014 140 2 | 04 015 015 140 2 | 04 015 016 140 2 |
| 1.75                                       | 04 175 009 140 2   | 04 175 010 140 2 | 04 175 012 140 2 | 04 175 014 140 2 | 04 175 015 140 2 | 04 175 016 140 2 |
| 2.0                                        | 04 020 009 140 2   | 04 020 010 140 2 | 04 020 012 140 2 | 04 020 014 140 2 | 04 020 015 140 2 | 04 020 016 140 2 |
| 2.25                                       | 04 225 009 140 2   | 04 225 010 140 2 | 04 225 012 140 2 | 04 225 014 140 2 | 04 225 015 140 2 | 04 225 016 140 2 |
| 2.5                                        | 04 025 009 140 2   | 04 025 010 140 2 | 04 025 012 140 2 | 04 025 014 140 2 | 04 025 015 140 2 | 04 025 016 140 2 |
| 2.75                                       | 04 275 009 140 2   | 04 275 010 140 2 | 04 275 012 140 2 | 04 275 014 140 2 | 04 275 015 140 2 | 04 275 016 140 2 |
| 3.0                                        | 04 030 009 140 2   | 04 030 010 140 2 | 04 030 012 140 2 | 04 030 014 140 2 | 04 030 015 140 2 | 04 030 016 140 2 |
| 3.25                                       | 04 325 009 140 2   | 04 325 010 140 2 | 04 325 012 140 2 | 04 325 014 140 2 | 04 325 015 140 2 | 04 325 016 140 2 |
| 3.5                                        | 04 035 009 140 2   | 04 035 010 140 2 | 04 035 012 140 2 | 04 035 014 140 2 | 04 035 015 140 2 | 04 035 016 140 2 |
| 3.75                                       | 04 375 009 140 2   | 04 375 010 140 2 | 04 375 012 140 2 | 04 375 014 140 2 | 04 375 015 140 2 | 04 375 016 140 2 |
| 4.0                                        | 04 040 009 140 2   | 04 040 010 140 2 | 04 040 012 140 2 | 04 040 014 140 2 | 04 040 015 140 2 | 04 040 016 140 2 |

| Shaft Length 1400mm, Guide Wire 0.014", RX |                    |                  |                  |                  |
|--------------------------------------------|--------------------|------------------|------------------|------------------|
| Balloon Diameter (mm)                      | Balloon Length(mm) |                  |                  |                  |
|                                            | 20                 | 24               | 25               | 30               |
| 1.5                                        | 04 015 020 140 2   | 04 015 024 140 2 | 04 015 025 140 2 | 04 015 030 140 2 |
| 1.75                                       | 04 175 020 140 2   | 04 175 024 140 2 | 04 175 025 140 2 | 04 175 030 140 2 |
| 2.0                                        | 04 020 020 140 2   | 04 020 024 140 2 | 04 020 025 140 2 | 04 020 030 140 2 |
| 2.25                                       | 04 225 020 140 2   | 04 225 024 140 2 | 04 225 025 140 2 | 04 225 030 140 2 |
| 2.5                                        | 04 025 020 140 2   | 04 025 024 140 2 | 04 025 025 140 2 | 04 025 030 140 2 |
| 2.75                                       | 04 275 020 140 2   | 04 275 024 140 2 | 04 275 025 140 2 | 04 275 030 140 2 |
| 3.0                                        | 04 030 020 140 2   | 04 030 024 140 2 | 04 030 025 140 2 | 04 030 030 140 2 |
| 3.25                                       | 04 325 020 140 2   | 04 325 024 140 2 | 04 325 025 140 2 | 04 325 030 140 2 |
| 3.5                                        | 04 035 020 140 2   | 04 035 024 140 2 | 04 035 025 140 2 | 04 035 030 140 2 |
| 3.75                                       | 04 375 020 140 2   | 04 375 024 140 2 | 04 375 025 140 2 | 04 375 030 140 2 |
| 4.0                                        | 04 040 020 140 2   | 04 040 024 140 2 | 04 040 025 140 2 | 04 040 030 140 2 |



### Princeton Medical Scientific Co., Ltd

Room 503, Building C, No.8 Fukang West Street,  
Yushu Industrial Park, Guangzhou Science City,  
Guangzhou, 510663, China

T: (86)20-82112462

F: (86)20-82119811

E: export-sales@pulinmed.com

www.pulinmed.com